

## Doctors with an iTero® scanner can send digital impressions directly to NDX Dental Art through iTero online.

## To send a file:

## iTero® Element™

- 1) Add **NDX Dental Art** (ID # **8056**) by calling iTero Support 800-577-8767.
- 2) Start a new Rx, and after entering patient name and case type, select NDX Dental Art from the Ship To menu.
- 3) If model-less, specify in the notes.
- 4) After scanning, press the Envelope icon, add signature, and press Confirm & Send.

## iTero® HD2.9 and Previous Models

- 1) Add NDX Dental Art (ID# 8056) by calling iTero Support at 800-577-8767 or from within your iTero online account.
- 2) Start a new Rx, enter required information and select NDX Dental Art from the Ship To menu.
- 3) If model-less, specify in the notes.
- 4) After completing the scan process, press the Send button.







iTero®

Should you have questions regarding the products, services and/or processes associated with submitting digital files to NDX Dental Art, please contact us directly at **800.444.3744** or **di-DentalArt@NationalDentex.com.** 

