## Carestream CS Scanner

Doctors with a Carestream CS Scanner can send digital impressions directly to NDX Keller through CS Connect or Direct E-mail.

## To send a file:

- 1) Complete the intraoral scan and open the patient's case in the DIS software. This software can operate in standalone mode or through Practiceworks.
- 2) Through the file button, in the upper right corner of the patient history window, select history open by clicking the book icon.
- 3) After highlighting the 3DIO folder, click the *Send to* button located toward the lower right corner of the window. The STL files will now be saved to the desktop.

## To submit to NDX Keller using CSCONNECT:

- a. After highlighting the appropriate 3DIO File from step 3, click the CS connect button. The web portal will open to login screen.
- b. Log into CS Connect.
- c. From the drop down menu, select NDX Keller.
  - Note: If this is the first case being sent to NDX Keller enter **Kellerstl-rc@kellerlab.com** in the email field. NDX Keller will appear in the drop-down menu moving forward.
- d. Attach files by clicking Browse and selecting the .dcm file from the patient's folder; then click Open.
- e. Add any pertinent information regarding the case to the comment box.
- f. Accept the terms and conditions by checking the box and click Submit.





## Carestream CS Scanner

Should you have questions regarding the products, services and/or processes associated with submitting digital files to NDX Keller, please contact us directly at **800.325.3056** or **Kellerstl-rc@kellerlab.com.**